

MESH Portable Household Atomizer NEBULIZER

JZ-492S



Description:

- 1) The fog is large, the atomized particles are fine, the particles can reach $3.7\mu\text{m}$, and the drug absorption is more complete;
- 2) Microporous spray sheet, atomization speed $\geq 0.25\text{ml}/\text{min}$;
- 3) Through piezoelectric components, low-noise design is more silent;
- 4) The minimum capacity of the medicine cup is about 0.5mL and the maximum is 10mL, and the residual liquid is less than or equal to 0.5ml. It is economical to use, and the micro-mesh cleaning function prevents the micro-mesh from clogging;
- 5) Two atomization methods can be selected according to different situations, namely, mask atomization (divided into human and child masks) and mouthpiece atomization;
- 6) The face mask is comfortable, skin-friendly and environmentally friendly;
- 7) Independent atomization components, easy to clean, one-click disassembly, easy to clean and disinfect, and avoid cross-infection;
- 8) When there is no liquid medicine in the cup or the voltage is insufficient, it will automatically shut down;
- 9) 2 kinds of power sources, convenient for home travel, use 2 AA alkaline batteries or use USB power adapter cable.

Technical Parameter

Mains connection	100-240V, 50-60HZ; 0.15A
Input	5V/1A;
Oscillation frequency	100kHz \pm 10%
Condensation rate	$\geq 0.25\text{ml}/\text{min}$.
Atomized particles	$\leq 5\mu\text{m}$
Medication flow rate approx	0.4ml/min.
Volume	max. 10ml
Product weight	115g \pm 10%

Packing Size

N.W.	99g (without battery)
G.W.	239g
Product size	67*42*116.5mm
Box size	140*110*60mm
QTY /G.W.	50pcs/12.7KG
CTN Size	574*300*310mm;
CBM	0.06m ³

Product List



CE RoHS FC