## ELECTRONIC

**Blood Pressure Monitor** 

JZ-256A



- 1) LCD liquid crystal digital display, the result shows high pressure / low pressure / pulse, large screen large characters, clear readings.
- 2) WHO sphygmomanometer classification indicates blood pressure health condition, convenient for users to accurately determine whether blood pressure is normal
- 3) Fully automatic intelligent pressurization and decompression, IHB arrhythmia detection
- 4) Year, month, day, hour, minute setting
- 5) 99 sets of measurement results memory for two people, can display the average reading of the ast 3 measurements, convenient for users to retrieve data comparison
- 6) Automatic shutdown after 1 minute of manual shutdown or no operation, easy and convenient operation
- 7) Blood pressure value unit Kpa and mmHg two display unit conversion (boot default unit is mmHg)
- 8) Overpressure protection function, automatic and rapid exhaustion when the air pressure exceeds 295mmHg-20ms, safe and
- 9) Error prompt function, convenient for users to find the cause of the error
- 10) Low-voltage detection function, low-voltage detection under any working condition,
- LCD display prompts low battery
- 11) measurable arm circumference: about 22~32cm
- 12) LCD display size: 41\*61mm











## **Technical Parameter**

recimiear raidificati	
Product name	Electronic blood pressure monitor (Arm type)
Display	LCD Digital Display
Memory array	99 sets of storage
Voice function	Optional
Measuring range	0mmHg~280mmHg
Maximum allowable	±3mmHg
error	
Operating Voltage	6V
Measurable arm	22cm~32cm
circumference	
Auto power off	Not Testing One Min Auto Shut Off

## **Packing Size**

N.W.	216g (without battery)
G.W.	398g
Product size	125*113*55mm
Box size	158*76*110mm
QTY/G.W.	30pcs/13KG
CTN Size	430*350*340mm;
СВМ	0.05m3

## **Product List**









